Surgical Silver Postoperative Hydrocolloid Dressings

McKESSON

McKesson surgical silver postoperative hydrocolloid dressings are sterile, single-use surgical dressings for the management of postoperative wounds healing by primary intent. Featuring an absorbent non-woven pad with ionic silver and a hydrocolloid adhesive structure, the dressing provides a barrier to bacterial penetration while supporting a moist wound environment. The flexible, waterproof, and breathable design allows patients to maintain daily activities while helping protect the incision site and support clinical wound management protocols.

Features and Benefits

- **ANTIMICROBIAL PROTECTION:** Contains ionic silver to provide a barrier against bacterial penetration, helping prevent colonization and proliferation of bacteria within the dressing for up to 7 days
- **BROAD-SPECTRUM EFFECTIVENESS:** Demonstrated *in vitro* effectiveness against the following Gram-positive bacteria: MRSA, MRSE, VRE; Gram-negative bacteria: *Escherichia Coli, Klebsiella Pneumoniae, Pseudomonas Aeruginosa;* Yeast: *Candida Albicans*
- **MOIST WOUND ENVIRONMENT:** Hydrocolloid layer forms a gel as fluid is absorbed, helping maintain a moist environment that supports wound healing
- **HIGHLY ABSORBENT PAD:** Absorbent non-woven central pad helps manage low to moderately exuding wounds while maintaining coverage over the incision
- **WATERPROOF & BREATHABLE BARRIER:** Flexible outer layer provides a waterproof yet breathable seal, allowing patients to shower while protecting the wound
- **SECURE, FLEXIBLE ADHESION:** Hydrocolloid adhesive layer conforms to wound contours and accommodates movement, supporting secure wear during patient activity
- **LOW-TRAUMA REMOVAL:** Designed to minimize trauma upon removal and allow intact removal without damaging sutures or closure systems
- **LEAKAGE CONTROL:** Windowed self-adhesive hydrocolloid contact layer helps minimize leakage while adhering to surrounding healthy skin
- **EXTENDED WEAR TIME:** Can remain in place for up to 7 days, depending on clinical assessment and exudate levels
- **NOT MADE WITH NATURAL RUBBER LATEX**





(16-915S: Front)




(16-915S: Back)

Surgical Silver Postoperative Hydrocolloid Dressings

McKESSON



WATERPROOF, BREATHABLE OUTER LAYER
Provides a waterproof barrier while remaining breathable to help protect the incision and allow patients to shower



(16-910S: Back)

ABSORBENT PAD WITH IONIC SILVER
Absorbs wound exudate while delivering ionic silver protection to help prevent bacterial colonization within the dressing for up to 7 days

SELF-ADHESIVE HYDROCOLLOID CONTACT LAYER
Adheres to surrounding skin to help secure the dressing, minimize leakage, and support low-trauma removal

GEL-FORMING HYDROCOLLOID TECHNOLOGY
Forms a gel as fluid is absorbed to help maintain a moist wound

Ordering Information

IMAGE	E1 #	MFR #	DESCRIPTION	SIZE	QTY
	1285224	16-910S	Surgical Silver Postoperative Dressing, Hydrocolloid, Single-use, Sterile	3.5 in x 4 in	10/bx 10 bxs/cs
	1285220	16-915S	Surgical Silver Postoperative Dressing, Hydrocolloid, Single-use, Sterile	3.5 in x 6 in	10/bx 10 bxs/cs
	1285221	16-925S	Surgical Silver Postoperative Dressing, Hydrocolloid, Single-use, Sterile	3.5 in x 10 in	10/bx 10 bxs/cs
	1285222	16-930S	Surgical Silver Postoperative Dressing, Hydrocolloid, Single-use, Sterile	3.5 in x 12 in	10/bx 10 bxs/cs
	1285223	16-935S	Surgical Silver Postoperative Dressing, Hydrocolloid, Single-use, Sterile	3.5 in x 14 in	10/bx 10 bxs/cs

Satisfaction Guaranteed
For complete details, please visit mms.mckesson.com/mckesson-brands.

Never Compromise on Care™

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

McKesson and the McKesson logo are registered trademarks of McKesson Corporation.

© 2026 McKesson Medical-Surgical Inc.

Rev. 00 06/26